



"The Clinical Advantage"™

ATOMLAB™

DOSE CALIBRATORS • WIPE TEST COUNTER

PROVEN
PERFORMANCE FOR
FAST, ACCURATE
MEASUREMENTS.



BRING SCIENCE,
TECHNOLOGY AND
APPLICATION
TOGETHER...WITH BIODEX

ISO 9001:2008
ISO 13485:2003
CERTIFIED

BIODEX
www.biodex.com
1-800-224-6339
Int'l 631-924-9000


Atomlab 500 and 500Plus **NEW** Features

- An enhanced “zero background” feature allows users to configure system background count time for 30, 60 or 100 seconds. This new feature changes the counting time when selecting the “zero background” button in all functions including the moly assay application.
- Expanded isotope libraries in both the dose calibrator and wipe test counter to include Lu-177 and Ra-223 isotopes. Libraries include updated dial values, isotope ROIs and efficiency settings.
- Half-Life Verification App included in the nuclear pharmacy menu. Simply select the isotope, choose two or three readings for the calculation, and set the time between readings. The results show the calculated half-life and the variance to the real half-life for the selected isotope.


Proven Accuracy & Quality Assurance




QA Dose Calibrator Screen – The extensive selection of quality assurance applications is readily available and simplifies regulatory requirements.



NEW Half-Life Verification Screen – Results show the calculated half-life and the variance to the published half-life for the selected isotope.




Constancy Data Screen – System stores and decay corrects multiple reference sources, compares the measured activity to the calculated activity and displays the percent variance. Data and source activity plots can be displayed, printed or stored.



Moly Assay Results Screen – Application simply guides you through the moly assay process making it quick and easy to perform, and displays the Moly concentration when test is complete. The system will prompt you if the results are not within regulatory limits.




Linearity Decay Plots Screen – Linearity tests can be performed using the traditional manual method or with the fully automated program featuring the capability to capture readings from a source, while automatically recording at specific intervals.




Dose Calculation Screen – All the information needed to draw doses efficiently is readily available on one screen. Easily performs pre- and post-decay calculations, volume calculations for specific times, and changes isotopes all with minimal screen touches. No calculators needed!

Quality Assurance Reports



Linearity Decay Plots
Report shows isotope, initial activity, elapsed time of readings, calculation method and variance, visually plotted.




Constancy Activity Plots
An Expanded Constancy test allows the user to measure a long-lived source and up to 14 routinely used isotope settings, all in a matter of seconds with just one touch of the screen. The system stores and decay corrects multiple reference sources and compares the measured activity to the calculated activity for the daily constancy test.

Full-size samples can be viewed at www.biodec.com/dosecalibrators

ATOMLAB™ 400 Dose Calibrator

...fast, accurate radionuclide activity measurements with performance that easily complies with the most stringent regulatory requirements.



Designed for facilities receiving unit doses including PET and BETA.

- Preprogrammed for 88 most commonly used radionuclides
- Automatic range selection; ranges up to 40 Curies of Tc-99m or 10 Curies of F-18
- Ultra-fast response
- Displays in Curies or Becquerels
- Remote Ionization Chamber
- RS-232 bi-directional serial communications port
- Self-Diagnostic Software
- Large, easy-to-read backlit LCD
- Small footprint economizes workspace
- Desktop or wall mount display
- Two-year warranty



◀ Atomlab™ 400 display can be mounted on a wall or placed on a desktop.

Atomlab™ 400

The Atomlab™ 400 provides fast, accurate radionuclide activity measurements with performance that easily complies with the most stringent regulatory requirements.

The unit is simple to operate. There is a routine list of ten pre-programmed isotopes plus another seven that are user selected from the library. The library contains 88 isotopes listed alphabetically, including Lu-177 and Ra-223. Four isotopes are displayed at a time.

Activity is displayed on a LCD panel in either Curie or Becquerel units. Background correction is performed at the touch of a button. Range selection is automatic.

Activity measurements are performed by a microprocessor-controlled electrometer located within the detector assembly of the ionization chamber. The chamber is shielded with .25" (6.3 mm) lead. It can be located up to eight feet away from the display unit. Chamber bias is generated within the display unit by an electronic high voltage supply, eliminating the need for expensive battery changes.

The RS-232 port enables the Atomlab 400 Dose Calibrator to communicate with most commercially available nuclear medicine management systems.

Specifications featured inside back cover.

BIODEX

www.biodex.com

1-800-224-6339

Int'l 631-924-9000

ATOMLAB™ 500 Dose Calibrator

Proven performance for fast, accurate measurements. One dose calibrator for all your requirements.

UPGRADE-ABILITY

The Atomlab™ 500 Dose Calibrator can be upgraded to include a fully functioning Atomlab™ Wipe Test Counter. The “Smart Display” recognizes the chamber / detector configuration and instantly reconfigures the screen to the appropriate icons.



Easy to use, touch-screen display with intuitive menus.

- Windows® Operating System
- Communicates with most commercially available NM management systems via Ethernet or Serial Port
- Pre-programmed for 89 most commonly used radionuclides; any 12 can be conveniently touch selected
- NEW Extensive Isotope Libraries include Lu-177 and Ra-223.
- Automatic range selection; ranges up to 100 Curies of Tc-99m or 25 Curies of F-18
- Ultra-fast response
- Displays in Curies or Becquerels
- USB ports accommodates a mouse, printing and software upgrades via memory devices.
- Report and label printers available
- Easy to use, large color touch screen display with intuitive menus
- Small footprint economizes workspace
- Upgradeable at any time to include a wipe test counter
- Desktop or wall mount display
- Two-year warranty



Standard Apps for Atomlab 500 include:

- Automated Quality Assurance Apps
 - Constancy and Expanded Constancy
 - Linearity and Auto Linearity
 - Accuracy
 - Geometry
- Nuclear Pharmacy Apps
 - Future dose and volume computation
 - Inventory control of 25 samples, correcting volume, activity and moly concentration
 - NEW Half-Life Verification
 - NEW Moly Assay

One dose calibrator that can be used for a wide variety of Nuclear Medicine, PET and radioimmunotherapy applications.

Atomlab™ 500

The Atomlab 500 can be used for a wide variety of nuclear medicine, PET and radioimmunotherapy applications, with proven performance for fast, accurate measurements. The system consists of a low pressure ionization chamber, electrometer with extraordinary linearity and an auto-ranging color touch screen display. Additionally, there are advanced, but easy-to-use programs for nuclear pharmacy, radiochemistry and radiochromatography.

Activity measurements are performed by the microprocessor controlled electrometer located within the chamber assembly. The chamber is shielded with .25" (6.3 mm) lead. It can be located up to eight feet away from the display unit. Chamber bias is generated by an electronic high voltage supply, eliminating the need for expensive battery changes.

Every element of the design and technical development will increase dose accuracy, department productivity and regulation compliance. The attractive and intuitive human interface guides the user through each operation. Software can easily be updated via the Biodex website or by using a convenient memory card. The touch-screen display can rest on a bench or mount on the wall of a hot lab, hot cell or laminar flow hood.


In addition to powerful self diagnostics, the Atomlab 500 includes an exclusive chamber monitoring technology to assure longer life and accuracy. Integrated pressure and temperature sensors feedback data so that the influence of gas pressure change will not effect an accurate reading.

OPERATION

The system is easy to use. There are 12 isotope selection touch keys pre-programmed for the most commonly used radionuclides. Any of those keys can be reprogrammed by the user for a desired isotope. There are 89 isotope-specific dial values listed in the library. Dial values can easily be changed if required.

Activity is displayed on the touch screen color display in either Curie or Becquerel units. Background correction is performed at the touch of a button. Range selection is automatic, from .01 microcurie to 100 Curies of Tc-99m or 25 Curies of F-18.

QUALITY ASSURANCE



The Atomlab 500 has been designed to make life easier. The extensive selection of quality assurance applications streamlines and simplifies hot lab administration requirements. The system stores and decay corrects multiple reference sources and compares the measured activity to the calculated activity for the daily constancy test.

Linearity tests can be performed in the traditional manual method or by a fully automated program that allows for readings from a source to be taken, and automatically recorded at specified intervals. The system will graph the results.

The attenuation tube test for linearity can be performed using software that will guide the user through the procedure, store the values and make all calculations.

Accuracy Test Dose Calibrator			
	Tested Source Co-57	Source ID BM06S-57-07-1	
	Current Activity	Calculated Standard	Acceptable Variance ± 2%
First	3.38 mCi	3.31 mCi	2.1%
Second	3.39 mCi	3.31 mCi	2.4%
Third	3.36 mCi	3.31 mCi	1.5%
Average	3.38 mCi	3.31 mCi	2.0%

ACCURACY

Atomlab Dose Calibrators have consistently proven to be highly accurate. Biodex and chamber manufacturer Sun Nuclear Corporation have participated in the isotope program sponsored by National Institute of Standards and Technology.

Each month a certified isotope is received from National Institute of Standards and Technology and is measured in the Atomlab™ Dose Calibrator, producing direct traceability to National Institute of Standards and Technology.

DATA DOWNLOAD

The Atomlab™ Data Manager is available as an option. The Windows™ based utility allows wipe test and dose calibrator QA results to be downloaded using a USB/serial converter. The results can be viewed and printed from the data manager software as required. In addition, information stored in the data manager can be exported into Microsoft® Excel or to department management systems.

BIODEX
www.biodex.com
1-800-224-6339
Int'l 631-924-9000

ATOMLAB™ 500 DOSE CALIBRATOR



086-341, Optional label printer for Atomlab 500 and Atomlab 500Plus

Sample Labels for Future Dose with Dose in Chamber

Future Dose Record

Patient: _____ Procedure: _____

Id #: _____ Kit: _____

Prepared By: _____ Lot #: _____

Isotope: Tc-99m

Calc Admin Dose: 106.3 mCi at 11:30 EST 06/01/2010

Measured Activity: 119. μ Ci at 10:33 06/01/2010

Administered By: _____

Radiopharmaceutical: Tc-99m

Calc Admin Dose: 106.3 mCi at 11:30 EST 06/01/2010

Measured Activity: 119. μ Ci at 10:33 06/01/2010

Radiopharmaceutical: Tc-99m

Calc Admin Dose: 106.3 mCi at 11:30 EST 06/01/2010

Measured Activity: 119. μ Ci at 10:33 06/01/2010

Save Expanded Constancy Data
Dose Calibrator

01/06/2014 Test Date

Biodex Model 086-332 Serial # 1008004p
Detector Serial # 5300111

Sealed Source	Source ID	Measured Activity	Calculated Standard	Acceptable Variance \pm 5%
Cs-137	BM06S-37	182.9 μ Ci	187.6 μ Ci	-2.5%

Isotope	Tc-99m	Cs-137	Ba-133	Co-57	Co-60	Cr-51	F-18
Current	297.	182.6	79.8	187.9	711.	53.4	3158.
Initial	298.	183.3	79.3	188.2	713.	53.5	3159.
Variance	0.6%	-0.3%	-0.6%	-0.3%	-0.3%	-0.2%	-0.3%

Staff Member CSS

Back Measure View Data Save Data

▲ Constancy Data

Moly Assay

Moly Assay is Complete

Make Note of the Concentration Before Exiting

Moly Breakthrough Reading 39.0 μ Ci

Tc-99m 3520. mCi

Moly Concentration 0.011 μ Ci Mo/mCi Tc-99m

01/07/2014 10:42

Back Measure

▲ Moly Assay

04/30/2008 Test Date

Geometry Test Data

Detector #1 Dose Calibrator

Biodex Model 086-330 Serial # 11111 Detector Serial # 1198883

Isotope Tc-99m Normalized Volume 3 ml

Volume	Measured Activity μ Ci	Calculated Activity μ Ci	Correction Factor (Not Applied) 0.95-1.05 Allowable
1 ml	246.	244.	0.99
2 ml	261.	244.	0.93
5 ml	244.	244.	1.00
10 ml	232.	243.	1.05
15 ml	262.	242.	0.92

Container Description

Type Syringe

Size 3 ml


Material Plastic

Technologist tom tom


Back Select Detector Select Isotope Delete Data Earlier Date Later Date Print

▲ Geometry Test

Easy to use large color touch screen with intuitive menus.



The Atomlab™ 500 display can be mounted on a wall or placed on a desktop



COMMERCIAL NUCLEAR PHARMACY

The Atomlab 500 Dose Calibrator features a nuclear pharmacy "Dose Calculation Screen" to meet the needs of a commercial nuclear pharmacy. All the information needed to draw doses efficiently is readily available on one screen. The feature easily performs pre- and post-decay calculations, volume calculations for specific times and isotope changes all with

minimal screen touches. No calculators needed!

Atomlab Dose Calibrators can read up to 100 Curies of Tc-99m, thereby eliminating the need for an aliquot preparation, a significant time saver every time a generator is milked.

INSTITUTIONAL NUCLEAR PHARMACY

The Atomlab 500 provides inventory control for 25 samples, storing and correcting the volume, activity, and moly concentration. The system will perform both volume and future dose calculations. In addition to inventory management, the Atomlab 500 provides quality assurance and record keeping functions. The inkjet printer allows hard copy records to be produced for all functions. The label printer allows the user to print labels for the syringe or vial.

RADIOCHROMATOGRAPHY

The radiopharmaceutical quality control program is exceptional. The Atomlab 500 performs all counting and calculations for paper chromatography tests, computing the percentages of free pertechnetate, hydrolyzed reduced Tc-99m and labeled radiopharmaceuticals.

RADIOCHEMISTRY

Up to seven ionization chambers, or six ionization chambers and a well counter can be connected via daisy chain to a single display. The activity in each detector can be selected and viewed from the single display.


COMMUNICATIONS

RS-232 port and two USB ports to communicate in real time with the most commercially available nuclear medicine management systems, connect to external monitor or upload software upgrades.

Specifications featured inside back cover.

EASY UPDATE


Your Atomlab 500 Dose Calibrator is upgradable. You can easily install software updates via the Biodex website or by using a convenient memory card.



ATOMLAB™ Wipe Test Counter

UPGRADE-ABILITY

The Atomlab™ Wipe Test Counter can be upgraded to include a fully functioning Atomlab™ Dose Calibrator simply by adding a new chamber. The “Smart Display” detects which chambers are connected and instantly reconfigures the screen to the appropriate icons.



Eliminate the tedium of wipe testing.

An extensive selection of quality assurance applications streamlines and simplifies hot lab administration requirements.

- Windows® Operating System
- Communicates with most commercially available NM management systems.
- USB Ports allow connection with external monitor and facilitate software upgrades
- 64 Channel MCA
- Adjustable wide window and individual isotopes
- 2x2 NaI drilled-well detector
- Remote shielded well
- Energy spectrums with individual ROI
- Ability to help identify isotopes causing contamination
- User-specific wipe locations and trigger levels
- Wipes that exceed trigger levels are immediately recognized
- Detailed wipe reports including cpm and dpm
- Wipe testing results stored
- Upgradable at any time to a dose calibrator by adding an ionization chamber
- Report and label printers available
- Two-year warranty



Standard Apps for Atomlab Wipe Test Counter include:

- Automated Quality Assurance Apps
 - Full Width at Half Max (FWHM)
 - Chi Square
 - Minimum Detectable Activity (MDA)

The Atomlab Wipe Test Counter is easy to use, easy to understand, fast and dependable. A color touch-screen display utilizing intuitive software eliminates the tedium of wipe testing. Simply perform a daily calibration and background count, then count the wipe for each predetermined location. Trigger levels can be set for any isotope at any location including 200 dpm for iodine. In seconds the system will determine if the location is above or below the user defined trigger level.

When performing a wipe test, the full spectrum is displayed. A wide window that includes the isotope energies expected in a particular department is set by the user. The efficiencies of the isotopes selected for the window can be either factory defaults or user determined using an integrated detector efficiency program. Individual isotope ROIs along with the wide window can be set. This feature helps identify the isotope(s) causing contamination.

Easy to use, easy to understand.

Up to 50 wipe locations can be entered as a restricted area, unrestricted area, sealed source or package. The results are displayed in dpm, cpm, μCi or Bq.

The system consists of a lead shielded 2" x 2" sodium iodide (NaI) well detector and a 64 channel multi-channel analyzer. The displayed energy range (spectrum) is 0-800 KeV, which is typically found in nuclear medicine departments.


The wipe counter is designed to meet or exceed all NRC (10 CFR 35.70, 10 CFR 20.1906 and 10 CFR 35.2067) and state wipe test requirements. There are automated programs for the quality assurance functions: calibration, FWHM, chi-square and minimum detectable activity (MDA).

Wipe test results and QA test data can be stored in memory and printed at any time.


Specifications featured inside back cover.

The Atomlab™ Wipe Test Counter display can be mounted on a wall or placed on a desktop.

Chamber shown with optional lead shield.



▲ Wipe Test Passed




▲ Wipe Test Failed



▲ QA Display (Wipe)



▲ MDA (Wipe)



▲ Chi Square

BIODEX

www.biodex.com
1-800-224-6339

Int'l 631-924-9000

ATOMLAB™ 500Plus Dose Calibrator

Dose Calibrator + Wipe Test Counter



Atomlab™ 500Plus...
comprehensive, easy-to-use
and feature rich.

One solution for all your molecular imaging needs.

The Atomlab 500Plus Dose Calibrator brings it all together - science, technology and application.

The Atomlab 500Plus combines the industry gold standard Atomlab 500 Dose Calibrator and Wipe Test Counter, offering you a complete and cost-effective solution for all of your molecular imaging needs. The intuitive Atomlab 500Plus provides fast, accurate radionuclide activity measurements with performance that complies with the most stringent regulatory requirements. It's comprehensive, easy-to-use and feature rich. The software mirrors the way you think and work. It guides when necessary, but does not

burden the advanced user. The touch-screen and easy-to follow on-screen prompts mean you won't need "cheat sheets."

When required, you can be confident that every operation is captured and documented. That documentation makes compliance a breeze.

The Atomlab 500 Dose Calibrator is upgradable. You can easily install software updates via the Biodex website or by using a convenient memory card.

- Windows® Operating System
- Communicates with most commercially available NM management systems via Ethernet or Serial Port
- Preprogrammed for 89 most commonly used radionuclides; any 12 can be conveniently touch selected
- Automatic range selection; up to 100 Curies of Tc-99m or 25 Curies of F-18
- Remote shielded ionization chamber
- Ultra-fast response
- Displays in Curies or Becquerels
- 64 Channel MCA
- Remote shielded well
- Adjustable wide window and individual isotopes
- Energy spectrums with individual ROI
- Ability to help identify isotopes causing contamination
- 2 x 2 NaI drilled-well detector
- User-specific wipe locations and trigger levels
- Wipes that exceed trigger levels are immediately recognized
- Detailed wipe reports including cpm and dpm
- Wipe testing results stored
- Intuitive software with extensive functionality

- USB port accommodates a mouse, printing and software upgrades via memory devices.
- Report and label printers available
- Easy-to-use, large, color touch-screen display with intuitive menus
- Small footprint economizes workspace
- Desktop or wall mount display
- Two-year warranty



Standard Apps for Atomlab 500Plus include:

- Automated Quality Assurance Apps
 - Constancy and Expanded Constancy
 - Linearity and Auto Linearity
 - Accuracy
 - Geometry
 - Full Width at Half Max (FWHM)
 - Chi Square
 - Minimum Detectable Activity (MDA)
- Nuclear Pharmacy Apps
 - Future dose and volume computation
 - Inventory control of 25 samples, correcting volume, activity and moly concentration
 - **NEW** Half-Life Verification
 - **NEW** Moly Assay

Specifications featured inside back cover.



▲ Atomlab 500Plus with additional chambers, wipe test counter and printer.

BIODEX
www.biodex.com
1-800-224-6339
Int'l 631-924-9000

ATOMLAB™ **500Plus** DOSE CALIBRATOR


ONLINE DEMONSTRATIONS & TRAINING AT YOUR DESK

Whether considering or learning, Biodex offers personal, online appointments for one-on-one sessions at **NO CHARGE**.


Experience the same level of detail as if we were on-site.

Currently available to customers in the USA

DOSE CALIBRATORS • THYROID UPTAKE • PULMONEX



▲ Wipe Test Passed



▲ Wipe Test Failed



▲ QA Display (Dose)



▲ Constancy Activity (Dose)

A complete Hot Lab Management System.

ACCESSORIES

MOLY ASSAY SHIELD




Moly Assay Shields are a convenient method to a unit dose. Whether working with syringes or vials, either unit will fit into the well chamber of any Atomlab Dose Calibrator.

086-435 Moly Assay Shield, Syringe, .3" lead

086-423 Moly Assay Shield, Vial, .3" lead

LINEATOR



The Lineator is a simple device to accurately and reliably verify the linearity of your dose calibrator. Test results are available in minutes, without waiting days for decay, making it feasible to perform a linearity test more often. Early identification can prevent problems before they occur.

086-509 Lineator

DOSE CALIBRATOR SHIELDING RINGS



The Dose Calibrator Shielding Rings offer an additional 2.25" of lead shielding around the remote chamber for working with 511 keV radionuclides, such as FDG F-18.

086-338 Shielding Rings, Interlocking, 2.25" lead (Fits Atomlab Dose Calibrators 400, 500 and 500Plus.) For additional protection from high energy activity

VIAL/SYRINGE DIPPER



This rugged, Vial/Syringe Dipper has a comfortable handle and it will hold 1 cc to 10 cc syringes or up to a 30 ml vial. Included with every Atomlab Dose Calibrator, the Vial/Syringe Dipper will also fit any well chamber with 2.5" x 10" interior dimension. The material used will not cause attenuation problems and is very resistant to breaking.

086-242 Dipper, Vial/Syringe


COPPER SYRINGE DIPPER



Designed for use with I-123 and In-111, the Copper Dipper removes variation in readings caused by attenuation differences from different materials and thicknesses used in syringes and vials.

086-242 Dipper, Vial/Syringe

The Atomlab™ 500Plus display can be mounted on a wall or placed on a desktop



▲ Linearity Decay (Dose)

BIODEX

www.biodex.com

1-800-224-6339

Int'l 631-924-9000



Atomlab™ 400

SPECIFICATIONS:

Isotope Selection Keys: Ten pre-programmed – Tc-99m, Tl-201, Co-57, Cs-137, I-131, In-111, Ga-67, Xe-133, I-123, and Mo-99; seven additional keys for user-set isotopes; two new isotope keys and a full alphabetical list of 88 isotopes.

Activity Range: 0.01 μ Ci to 40 Ci (.0004 MBq to 1500 GBq) of Tc-99m or 10 Ci of F-18.

Energy Range: 25 keV to 3 MeV photons

Response Time: One to two seconds for doses greater than 200 uCi; three seconds for doses greater than 20 uCi; 50-100 seconds below 20 uCi of Tc-99m with default threshold, threshold adjustable to reduce counting time

Detector Linearity: $\pm 1\%$ or 0.2 μ Ci, whichever is greater

Electrometer Linearity: $\pm 1\%$ or 0.2 μ Ci, whichever is greater

Electrometer Accuracy: $\pm 1\%$ or 0.2 μ Ci, whichever is greater

Overall Accuracy: $\pm 3\%$ or 0.3 μ Ci, whichever is greater; overall accuracy is affected by such factors as the accuracy of the specific source calibration, geometric variations due to sample volume or configuration, detector linearity, electrometer accuracy and readout accuracy

Repeatability: $\pm 0.3\%$ above 1 mCi short term (24 hr); 1% long term (one yr); exclusive of background

Detector: Well-type pressurized ionization chamber, with Argon fill gas; well opening 2.75" (7 cm), well depth 10.25" (26 cm)

Detector Shielding: .25" (6.3 mm) lead on all sides except top well opening; supplementary shielding available

Power Requirements: 100 to 240 VAC, 0.6 – 0.3 amps, auto switching

Display Unit:

Dimensions: 6.75" w x 6" depth x 5" h (17.1 x 15.3 x 12.7 cm)

Weight: 3.6 lb (1.64 kg); desktop or wall mountable

Detector Unit:

Dimensions: 6" dia x 15.5" h (15.24 x 39.37 cm)

Well I.D.: 2.75" dia x 10.5" h (7 x 26.7 cm)

Well I.D. with Liner: 2.5" dia x 10.25" h (6.35 x 26 cm)

Lead Shielding: .25" thick (6.3 mm)

Weight: 35 lb (16 kg)

Certification: ETL listed to UL 60601-1 and to CAN/CSA C22.2 No. 601-1-M90

IEC 60601-1, IEC 60601-1-4 and IEC 60601-1-2 and CE marked

Warranty: Two years parts and labor



086-335 Dose Calibrator, Atomlab™ 400, 100-240 VAC

Includes: Smart Display, ionization chamber, RS-232 port, vial/syringe dipper and well insert.



Atomlab™ 500

SPECIFICATIONS:

Display: LCD Touch Panel 6.5" x 5", function keys are displayed for the operation being performed

Dimensions: 9.5" w x 12" depth x 12" h (24.1 x 30.5 x 30.5 cm)

Weight: 6.3 lb (2.9 kg); desktop or wall mountable

Line Voltage: 100 to 240 VAC, auto selectable by the power supply

Line Frequency: 50/60 Hz, detachable line cord, built-in EMI filter and transient suppression

Auxiliary Port: Two USB ports, one RS-232 port

Isotope Selection Keys: Twelve pre-programmed – Tc-99m, Tl-201, I-123, I-131, Cs-137, Co-57, Xe-133, Ga-67, In-111, F-18, Y-90s, Ba-133, 25 user-defined isotopes and a full alphabetical list of 89 isotopes.

Activity Range: 0.01 uCi to 100 Ci (.0004 MBq to 3700 GBq) of Tc-99m or 25 Ci of F-18

Energy Range: 25 keV to 3 MeV photons

Response Time: One to two seconds for doses greater than 200 uCi; three seconds for doses greater than 20 uCi; 50-100 seconds below 20 uCi of Tc-99m with default threshold; threshold adjustable to reduce counting time

Detector Linearity: $\pm 1\%$ or 0.2 μ Ci, whichever is greater

Electrometer Linearity: $\pm 1\%$ or 0.2 μ Ci, whichever is greater, up to 40 curies of Tc-99m, $\pm 1.5\%$ up to 100 curies of Tc-99m

Electrometer Accuracy: $\pm 1\%$ or 0.2 μ Ci, whichever is greater

Overall Accuracy: $\pm 3\%$ or 0.3 μ Ci, whichever is greater; overall accuracy is affected by such factors as the accuracy of the specific source calibration, geometric variations due to sample volume or configuration, detector linearity, electrometer accuracy and readout accuracy

Repeatability: $\pm 0.3\%$ above 1 mCi short term (24 hr); 1% long term (one yr)

Detector: Well-type pressurized ionization chamber, with Argon fill gas; well opening 2.75" (7 cm), well depth 10.25" (26 cm). Up to seven chambers can be serially connected to one display

Detector Shielding: .25" (6.3 mm) lead on all sides except top well opening; supplementary shielding available

Power Requirements: 100 to 240 VAC, 0.38 – 0.15 amps, auto switching

Detector Unit:

Dimensions: 6" dia x 15.5" h (15.24 x 39.37 cm)

Well I.D.: 2.75" dia x 10.5" h (7 x 26.7 cm)

Well I.D. with Liner: 2.5" dia x 10.25" h (6.35 x 26 cm)

Lead Shielding: .25" lead (6.3 mm)

Weight: 35 lb (16 kg)

Certification: ETL listed to UL 60601-1 and to CAN/CSA C22.2 No. 601-1-M90

IEC 60601-1, IEC 60601-1-4 and IEC 60601-1-2 and CE marked


Warranty: Two years parts and labor



086-330 Dose Calibrator, Atomlab™ 500, 100-240 VAC

Includes: Smart Display, ionization chamber, RS-232 port, vial/syringe dipper and well insert.

Atomlab Dose Calibrators utilize low-pressure chambers and are shipped air or ground as standard goods.



Atomlab™ WIPE TEST COUNTER

SPECIFICATIONS:

Display: LCD Touch Panel 6.5" x 5", function keys are displayed for the operation being performed

Dimensions: 9.5" w x 12" h x 12" d (24.1 x 30.5 x 30.5 cm)

Weight: 4.2 lb (1.9 kg)

Line Voltage: 100 to 240 VAC, auto selectable by the power supply

Line Frequency: 50/60 Hz, detachable line cord, built-in EMI filter and transient suppression

Auxiliary Port: Two USB ports, one RS-232 port

Memory: Stores wipe, calibration, background, high voltage, isotope specification, isotope efficiency, Chi-Square testing results, technologist list, wipe locations list and latest MDA calculation. Results can be displayed and printed.

Preset Radionuclides: 27 including Tc-99m, Co-57, Cs-137, Ga-67, Tl-201, I-123, I-125, I-131, In-111, F-18

WELL COUNTER

Dimensions: 6" d x 11" h (15.24 x 27.9 cm)

Weight: 29 lb (3.2 kg)

Detector: 2" x 2" NaI (TI) integral line scintillation detector with a 0.75" dia x 1.44" depth well (1.9 x 3.7 cm)

Style: Remote Detector

Channels: 64

MCA: Integral to Well Counter

Spectral Resolution: FWHM 10%

Count Rate: (Maximum) 30,000 cps

Lead Shielding: 0.5" (1.2 cm) integral lead shield

Certification: ETL listed to UL 60601-1 and to CAN/CSA C22.2 No. 601-1M90, IEC 60601-1, IEC 60601-1-4 and IEC 60601-1-2 and CE marked

Warranty: Two years parts and labor



086-331 Atomlab™ Wipe Test Counter
Includes: Smart Display, well counter and RS-232 port.



Atomlab™ 500Plus

SPECIFICATIONS:

For specifications and features, see the Atomlab™ 500 Dose Calibrator and the Atomlab Wipe Test Counter. The Atomlab™ 500Plus incorporates the features and specifications of both units in one compact footprint.



086-332 Dose Calibrator, Atomlab™ 500Plus, 100-240 VAC
Includes: Smart Display, ionization chamber, well counter, RS-232 port, vial/syringe dipper and well insert.

RELATED PRODUCTS

Related: Dose Calibrator

Data Manager

086-333 Software, Atomlab 500 Data Manager

086-336 Chamber, Dose Calibrator

075-594 Chamber, Wipe Test

086-338 Shielding Rings, Interlocking, 2.25" lead
For additional protection from high energy activity

086-509 Lineator

086-435 Moly Shield, Syringe, .3" lead

086-243 Copper Dipper, Vial/Syringe

086-423 Moly Shield, Vial, .3" lead

086-341 Printer, Label, Dymo
Includes one roll 086-343 labels; two rolls 086-344 label

086-339 Printer, Ink Jet (report)

086-334 Cable, European to Wall Outlet

Replacement:

086-242 Vial/Syringe Dipper

086-241 Well Insert

086-343 Label, Blank, Lg, 300/roll
(For Dymo Printer; 086-341)

086-344 Label, "Radioactive", Sm, 260/roll, 2/pkg
(For Dymo Printer; 086-341)

Related: Wipe Test Counter

063-139 Rod Source, Cs-137, Calibrated, 0.1 µCi

075-596 Lead Shield, Wipe Test Chamber

006-350 Wipe Test Kits, 500/pkg

086-342 Well Liners, Disposable, 100/pk

An industry exclusive two-year warranty is standard.

For physics tests, sample reports, testimonials and warranty information, visit us on the web www.biodex.com/dosecalibrators

BIODEX
www.biodex.com
1-800-224-6339
Int'l 631-924-9000

BIODEX

Biodex Medical Systems, Inc.
20 Ramsey Road, Shirley, New York, 11967-4704

CHANGE SERVICE REQUESTED

**SUBSCRIBE to
BIODEX**

- NEWSLETTERS
- BLOGS
- EMAIL

www.biodex.com/eregister



Follow Biodex-Nuclear Medicine
Let's get a discussion going...




FN: 15-149 4/15

"The Clinical Advantage"™

Atomlab 960 Thyroid Uptake System

A complete, mobile, self-contained Medical Spectrometer System

Unique-positioning LED for accurate thyroid centering,
a first in Thyroid Uptake System design.



Optional **DICOM** Compliant
Software Program For
Atomlab™ Thyroid Uptake
System

www.biodex.com/thyroiduptake

BIODEX