

AIR-CUSHIONED FACE MASK



The SoftSeal air-cushioned face mask has a pre-filled air-cushion (medium inflation pressure) that molds to the contour of the patient's face providing a tight, leak-proof and comfortable seal. Cushion pressure is easily adjusted by inserting a standard syringe without needle in the two-way valve located on the mask bottom. The flexibility of cushion pressure allows optimum surface contact for every patient. Transparency of the entire mask allows continuous visual identification of patient's vital signs.

SOFTSEAL™ Air-Cushioned Face Mask without Injection Port:
132-685 Face Mask without Injection Port, Adult ☺

AIR-CUSHIONED FACE MASK



With direct dose administration

In addition to the patient comfort provided by the soft, pliable air cushion of the disposable face mask, the direct dose feature reduces patient anxiety because no needle is used. The syringe luer locks tightly to the direct dose tube, away from the patient's face. The one-way valve prevents the gas from reversing through the mask luer connector, allowing the technologist to disconnect immediately after injecting.

SOFTSEAL™ Air-Cushioned Face Mask with Direct Dose Administration:
132-698 Face Mask with Direct Dose Administration, Adult ☺



SOFTSEAL™ Air-Cushioned Face Mask with Injection Port:
132-695 Face Mask with Injection Port, Adult ☺

FACE MASK HARNESS



This traditional Face Mask Harness is made of soft rubber and can be adjusted for any head size. The square rubber base is comfortably positioned on the back of the patient's head and the harness tails are brought around to the front. The small holes fit snugly on the hook ring and hold the mask firmly in place.

139-677 Face Mask Harness, Adult

800.409.2688



technologyimaging.com