

DISPOSABLE SHOE COVERS



Disposable, conductive shoe covers are made of non-skid Tyvek with attached conductive strip to ground static electricity. White color, fits all sizes.

119-200 Shoe Covers, Disposable, 10 pr

DISPOSABLE COVERALL



Disposable Coveralls were designed with comfort and safety in mind. They are made of Tyvek, a tough material that is virtually tear-free and puncture proof. The coveralls actually exceed the standards established by the Federal Flammable Fabrics Act.

The one-piece design, complete with hood and zipper closure, eliminates the problem of garment maintenance and decontamination.

118-125 Coverall, Disposable, Large

COLLIMATOR PROTECTORS



Shown with Wall Dispenser, 114-017

Disposable plastic Collimator Protectors are designed to keep the face of the scintillation camera collimators from becoming contaminated by infectious fluids.

Protector sheets are backed with an adhesive for quick and easy application. After use,

peel away the contaminated sheet and discard. Apply a new protector between each patient.

Protectors are packaged in perforated rolls. The optional dispenser is designed of heavy gauge polished stainless steel with an aluminum rod. Mount to any wall or door for easy access.

SPECIFICATIONS:

Dimensions: 17" x 17" (43 x 43 cm)

Weight: 9 lb (4 kg)

114-177 Collimator Protectors, 250/roll

Related:

114-017 Wall Dispenser, Collimator Protectors
Wall mounted

ABSORBENT PAPER



Highly absorbent with plastic lining for fast, safe clean up

Protect any work surface with plastic lined Absorbent Paper. Spills and splashes are immediately contained without messy clean up. Available in pre-cut sheets or by the roll, Absorbent Paper is efficient for any work surface.

Together with the optional paper dispenser, the 300-foot Absorbent Paper Roll can be conveniently located to any wall or bench. It works like a common paper towel dispenser, yet is extra-strong to support the oversized 11 lb roll.

SPECIFICATIONS:

033-300 Absorbent Paper Roll

Dimensions: 300' l x 20" w (91.5 m x 51 cm)

033-013 Absorbent Paper Sheets

Dimensions: 17" l x 13" w (43 x 33 cm)

033-300 Absorbent Paper, Roll, 300 ft

033-013 Absorbent Paper, Sheets, 50/pkg

Related:

046-275 Dispenser, Absorbent Paper Roll
Wall mounted

800.409.2688



Clear Solutions, Uncommon Service.

technologyimaging.com