

DOSE DRAWING SYRINGE SHIELDS

For drawing a dose from a shielded vial



Dose Drawing Syringe Shield (glass non-replaceable)



Dose Drawing Syringe Shield with replaceable glass

Dose Drawing Syringe Shields reduce hand exposure when drawing doses from vial shields. The barrel of the syringe is constructed of optically clear high density (5.6) lead glass that offers complete 360° visibility and reduces radiation exposure from Tc-99m by 98.4%. At the end of the lead glass barrel is a .125" lead flange that provides additional hand shielding between a vial shield and the syringe shield.

Dose Drawing Syringe Shields:

Accommodates luer and non-luer lock syringes

- 007-661** Syringe Shield, 3 cc
- 007-663** Syringe Shield, 5 cc and 6 cc
- 007-665** Syringe Shield, 10 cc and 12 cc

Dose Drawing Syringe Shields:

Accommodates B-D Safety-Lok and Monoject safety syringes

- 007-645** Syringe Shield, Safety Syringe, 3 cc
- 007-646** Syringe Shield, Safety Syringe, 5 cc
- 007-647** Syringe Shield, Safety Syringe, 10 cc

The shield features quick and smooth syringe insertion with an O-ring seal and anti-roll cap. Dose Drawing Syringe Shields for safety syringes facilitate two thumbscrews to secure the syringe.

The shields are constructed with an elastomer material that acts as a shock absorber between the metal and glass, improving the resistance to breakage.

Dose Drawing Syringe Shields with Replaceable Glass:

Accommodates luer and non-luer lock syringes

- 007-691** Syringe Shield, 3 cc
- 007-693** Syringe Shield, 5 cc and 6 cc
- 007-695** Syringe Shield, 10 cc and 12 cc

Replacement Glass for Dose Drawing Syringe Shields:

Accommodates luer and non-luer lock syringes

- 127-691** Syringe Shield Replacement Glass, 3 cc
- 127-693** Syringe Shield Replacement Glass, 5 cc and 6 cc
- 127-695** Syringe Shield Replacement Glass, 10 cc and 12 cc

800.409.2688



technologyimaging.com