

# Flood Source Shielding “Radshield”



Radshield™ was designed to provide an easier alternative for transferring your source to the camera with minimal exposure without moving a 60 plus pound case around your Nuclear Medicine Department. While weighing only 28 pounds it provides similar radiation dose reduction as the typical hard case, by shielding the surface of the source not an oversize case. An optional integrated hard case or Radscooter is available from RadQual if desired.

	<u>Only Radshield™</u>	<u>Competitor's Hard Case</u>
		(millirem/hour/millicurie)
@ contact	0.1	0.07
@ 30 cm	0.04	0.03
@ 1 m	0.01	0.01



**tis**

Clear Solutions. Uncommon Service.

technologyimaging.com  
PO Box 3589  
Youngstown, OH 44513  
330.758.7800  
800.409.2688

Based on a 15 mCi Co-57 source using a Bicron MicroRem meter with Co-58/Co-58 levels below 0.08% combined.