

SHARPS CONTAINER SHIELDS

Key-lock design meets OSHA standards



Transfer the full Sharps Container to the Decay Module for syringe decay storage.

Sharps Container Shields are a simple, safe and convenient way to dispose of used syringes that may contain low-energy gamma radiation residue. The shields are constructed of steel lined with .125" (.32 cm) lead. They feature a top with sliding pott and a swing handle that allows easy transport.

Two models are available for one sharps container:

- Model 039-325 accommodates small and medium Monoject Sharps Containers.
- Model 039-335 accommodates Sharpstainer* and B-D Guardian Sharps Containers.

One model is available for two sharps containers:

- Model 039-350 accommodates two small Monoject Sharps Containers simultaneously for decay rotation.

The sharps containers are stacked inside the shield. When the top container is filled, the bottom decayed container is removed for proper waste disposal. The top container is then shifted to the lower position for decay while a fresh container is placed on top for immediate use. Spring clamps hold containers securely in place.



SPECIFICATIONS:

039-325 Sharps Container Shield
Shield uses 039-338 and 039-387
Monoject Sharps containers

Dimensions: 7.5" l x 11.75" w x 10" h
(19 x 30 x 25 cm)

I.D.: 7" l x 11.25" w x 9.5" h
(18 x 29 x 24 cm)

Lead Shielding: .125" thick (.32 cm)

Security: Key-locked

Finish: Powder coat

Weight: 48 lb (22 kg)



039-335 Sharps Container Shield
Shield uses 039-341 and
039-321 Sharps containers

Dimensions: 7.5" l x 12.5" w x 13" h
(19 x 32 x 33 cm)

I.D.: 6.875" l x 10.75" w x 12.5" h
(17.5 x 27 x 32 cm)

Lead Shielding: .125" thick (.32 cm)

Security: Key-locked

Finish: Powder coat

Weight: 53 lb (24 kg)



039-350 Sharps Container Shield, Dual
Shield uses two 039-338 Monoject Sharps containers

Dimensions: 7.375" l x 12" w x 16.75" h
(19 x 30.5 x 42.5 cm)

I.D.: 7" l x 11.5" w x 14.875" h
(18 x 29 x 38 cm)

Lead Shielding: .125" thick (.32 cm)

Finish: Powder coat

Weight: 60 lb (27 kg)



Sharps Container Shields:

- 039-325** Sharps Container Shield, .125" lead
For Monoject containers 039-338 and 039-387
- 039-335** Sharps Container Shield, .125" lead
For Sharpstainer* container 039-341 and
B-D Guardian container 039-321
- 039-350** Sharps Container Shield, Dual, .125" lead
For Monoject containers 039-338

Related:

- 039-338** Monoject Container, 4 qt., 10/pkg
Fits 039-325, 039-326, 039-350 and 039-330
- 039-387** Monoject Container, 8 qt., 10/pkg
Fits 039-325, 039-326 and 039-330
- 039-341** Sharpstainer* Container, 6.2 qt., (#182),
12/pkg
Fits 039-335
- 039-321** B-D Guardian Container, 8 qt., 8/pkg
Fits 039-335

800.409.2688



Clear Solutions. Uncommon Service.

technologyimaging.com