

DISPOSABLE MOUTHPIECE



This mouthpiece has a soft flange, formed tooth-grip protrusions, and fits snugly between the patient's lips and teeth to form a tight seal.

130-550 Mouthpiece, Disposable ©

TRU-FIT DISPOSABLE MOUTHPIECE



Contoured to the natural shape of the mouth, this mouthpiece can be held gently, but securely — eliminating the need to clamp down or the possibility of “popping” out of the patient's mouth and interrupting the study. The soft plastic material is transparent so that any obstruction can be immediately detected by the technologist. Pleasantly scented

with vanilla, this mouthpiece makes the procedure a little more comfortable for the patient.

130-551 Mouthpiece, Disposable, Vanilla Scented ©

DISPOSABLE NOSE CLIP



Economical, this nose clamp is all plastic for single patient use. Used with any mouthpiece.

130-100 Nose Clip, Disposable, 100/pkg ©

DISPOSABLE BACTERIA FILTER



Bacteria filters are used to reduce the possibility of cross contamination. The single-use filter is placed in line between a delivery tube and disposable mouthpiece or face mask. Electrostatically charged filter media is 99.9% effective in bacteria/virus retention while maintaining low breathing resistance. Two filters can be piggy-backed together for potential high-risk studies.

132-750 Bacteria Filter ©

800.409.2688



technologyimaging.com