

























COLOR-CODED SYRINGE LABELS



Labels quickly identify the radiopharmaceutical

 MDP 212-001 <i>Purple</i>	 TI 212-009 <i>Reflex Blue</i>	 CARDIOLITE 212-019 <i>Pink</i>
 MAA 212-002 <i>Light Orange</i>	 212-011 <i>Unmarked White</i>	 CHOLETEC 212-023 <i>Orange</i>
 SC 212-003 <i>Yellow</i>	 HIDA 212-012 <i>Brown</i>	 Tc RBC 212-026 <i>White</i>
 DTPA 212-004 <i>Olive Green</i>	 HDP 212-013 <i>Purple</i>	 MYOVUE 212-027 <i>White</i>
 GH 212-005 <i>Red</i>	 MIBI 212-014 <i>White</i>	 In-111 212-030 <i>White</i>
 Tc04 212-006 <i>Pantone Purple</i>	 DMSA 212-015 <i>Green</i>	 Saline 212-031 <i>White</i>
 PYP 212-007 <i>Blue</i>	 MAG-3 212-016 <i>Blue</i>	 FDG 212-032 <i>White</i>
 Ga 212-008 <i>Orange</i>	 I-131 212-017 <i>White</i>	
 NeuroLite 212-028 <i>White</i>		

This color-coded system has been designed to quickly identify radiopharmaceuticals. The .5" (1.3 cm) diameter circles are colored, labeled with the radiopharmaceutical name, and easily affixed to the plunger top of syringes. The colors have been selected as a standard most representative of radiopharmaceutical manufacturers' packaging schemes.

Color-Coded Syringe Labels (500/roll):

212-001	Label, MDP
212-002	Label, MAA
212-003	Label, SC
212-004	Label, DTPA
212-005	Label, GH
212-006	Label, Tc04
212-007	Label, PYP
212-008	Label, Ga
212-009	Label, TI
212-011	Label, (blank)
212-012	Label, HIDA
212-013	Label, HDP
212-014	Label, MIBI
212-015	Label, DMSA
212-016	Label, Mag-3
212-017	Label, I-131
212-019	Label, Cardiolite
212-023	Label, Choletec
212-026	Label, Tc RBC
212-027	Label, Myoview
212-028	Label, NeuroLite
212-030	Label, In-111
212-031	Label, Saline
212-032	Label, FDG