

SHIELDED SYRINGE HOLDER

Accommodates shielded and unshielded syringes



This Shielded Syringe Holder will accommodate unshielded syringes and syringes in a syringe shield.

The Syringe Holder is constructed of lead shielding, encased in steel. The shielding tapers from .25" - .5" lead. The large diameter base ensures stability.

For added protection, store loaded syringes in the Shielded Syringe Holder

SPECIFICATIONS:

Dimensions: 6.5" h (16.5 cm)

I.D.: .84" dia. x 5.6 h (2.1 x 14.3 cm)

Lead Shielding: .25" - .5" thick (.64 cm - 1.3 cm)

Accommodates Syringe Shields:

Pro-Tic II: 1 cc, 3 cc, 5 cc

Pro-Tic III: 1 cc, 3 cc, 5 cc

Pro-Tic IV: 1 cc

Weight: 6 lb (2.7 kg)

009-205 Syringe Holder, Shielded

UNIT DOSE PIG WALL RACK



Improve lab safety, efficiency and organization with the new Unit Dose Pig Wall Rack.

The wall rack improves work space with its pigeonhole design and reduces unnecessary handling. Unit doses can be identified at a glance. The rack's sturdy construction will hold up to 25 unit dose pigs.

Holds pigs up to 2.25" diameter.


SPECIFICATIONS:

Dimensions: 17" w x 18" h (43.2 x 45.7 cm)

008-400 Wall Rack, Unit Dose Pig

SYRINGE SHIELD HOLDER

Protect your investment



Don't let your syringe shields roll around on the counter. The Syringe Shield Holder offers a means of protecting syringe shields from scratches or misplacement while freeing up extra work space. The Syringe Shield Holder will support up to eight shields

and is counterbalanced to prevent tipping. Know exactly where syringe shields are when you need them.

SPECIFICATIONS:

Dimensions: 7.5" w x 4" depth x 6" h (19 x 10.2 x 15.2 cm)

Weight: 5 lb (2.3 kg)

007-999 Syringe Shield Holder

SYRINGE RECAPPER

Syringe recapping device



Don't put yourself at risk with an accidental needle stick! The Syringe Recapper is a safe and inexpensive way to protect yourself when recapping a used syringe. Used either hand-held or placed on a flat surface, such as a procedure tray, the Recapper is made of a lightweight plastic that is easily carried anywhere.

008-300 Syringe Recapper

3.25" l x 2.75" w (8.3 x 7 cm)

SHIELDED SYRINGE CARRIERS

Two sizes, each offered in .125" and .25" lead shielding for added protection



Shielded Syringe Carriers, small and large

Shielded Syringe Carriers reduce exposure while storing or transporting radioactive material. The overlapping lid design with snap-latch closure prevents streaming. There are two sizes and thicknesses of lead to choose from. The ends of the carriers are double thick to reduce the exposure from the ends of syringes.

SPECIFICATIONS:

001-179 Shielded Syringe Carrier, Large

Dimensions:
 I.D.: 8.25" l x 3" w x 2.9" h (21 x 7.6 x 7.4 cm)
 O.D.: 9.25" l x 4.4" w x 3.5" h (24 x 11.2 x 8.9 cm)
Lead Shielding:
 Sides, top and bottom: .125" thick (.32 cm)
 Ends: .25" thick (.64 cm)
 Weight: 11.3 lb (5.1 kg)

001-181 Shielded Syringe Carrier, Small

Dimensions:
 I.D.: 8" l x 1.9" w x 1.97" h (20.5 x 4.8 x 5 cm)
 O.D.: 9.25" l x 3.4" w x 2.6" h (23.5 x 8.6 x 6.6 cm)
Lead Shielding:
 Sides, top and bottom: .125" thick (.32 cm)
 Ends: .25" thick (.64 cm)
 Weight: 7.5 lb (3.4 kg)

001-182 Shielded Syringe Carrier, Small

Dimensions:
 I.D.: 7.5" l x 1.7" w x 1.7" h (19 x 4.3 x 4.3 cm)
 O.D.: 9.25" l x 3.4" w x 2.6" h (23.5 x 8.6 x 6.6 cm)
Lead Shielding:
 Sides, top and bottom: .25" thick (.64 cm)
 Ends: .5" thick (1.3 cm)
 Weight: 11 lb (4.9 kg)

001-180 Shielded Syringe Carrier, Large

Dimensions:
 I.D.: 7.7" l x 2.6" w x 2.7" h (19.6 x 6.6 x 6.9 cm)
 O.D.: 9.5" l x 4.4" w x 3.5" h (24 x 11.2 x 8.9 cm)
Lead Shielding:
 Sides, top and bottom: .25" thick (.64 cm)
 Ends: .5" thick (1.3 cm)
 Weight: 17 lb (7.7 kg)

Shielded Syringe Carrier, .125" lead:

- 001-181** Syringe Carrier, Small
- 001-179** Syringe Carrier, Large

Shielded Syringe Carrier, .25" lead:

- 001-182** Syringe Carrier, Small
- 001-180** Syringe Carrier, Large